

Phospho-PDGFR- β (Y751) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15725	SYNONYMS
SIZE	100 ug	PDGFRB; PDGFR; PDGFR1; Platelet-derived growth factor receptor beta; PDGF-R-beta; PDGFR-beta; Beta platelet-derived growth factor receptor; Beta-type platelet-derived growth factor receptor; CD140 antigen-like family member B; Platelet-deri
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	PDGFRB
PROTEIN NAME	Platelet-derived growth factor receptor beta
IMMUNOGEN	Synthesized peptide derived from human PDGFR- β around the phosphorylation site of Y751.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	IHC,IF,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.