

## PDGFR- $\beta$ Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX15727	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	PDGFRB; PDGFR; PDGFR1; Platelet-derived growth factor receptor beta; PDGF-R-beta; PDGFR-beta; Beta platelet-derived growth factor receptor; Beta-type platelet-derived growth factor receptor; CD140 antigen-like family member B; Platelet-deri
<b>HOST</b>	Rabbit	

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	PDGFRB
<b>PROTEIN NAME</b>	Platelet-derived growth factor receptor beta
<b>IMMUNOGEN</b>	Synthesized peptide derived from human PDGFR- $\beta$ around the non-phosphorylation site of Y1021.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human,Mouse,Rat
<b>TESTED APPLICATIONS</b>	WB,IHC,IF,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**