

ATP5G1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15729	SYNONYMS
SIZE	100 ug	ATP5G1; ATP synthase lipid-binding protein;
HOST	Rabbit	mitochondrial; ATP synthase proteolipid P1; ATPase protein 9; ATPase subunit c

PRODUCT INFORMATION

GENE NAME	ATP5G1
PROTEIN NAME	ATP synthase lipid-binding protein mitochondrial
IMMUNOGEN	Synthesized peptide derived from the Internal region of human ATP5G1.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.