

ERGIC-3 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15730	SYNONYMS
SIZE	100 ug	ERGIC3; C20orf47; ERV46; SDBCAG84; CGI-54;
HOST	Rabbit	Endoplasmic reticulum-Golgi intermediate compartment protein 3; Serologically defined breast cancer antigen NY-BR-84

PRODUCT INFORMATION

GENE NAME	ERGIC3
PROTEIN NAME	Endoplasmic reticulum-Golgi intermediate compartment protein 3
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human ERGIC-3. at AA range: 280-360

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse
TESTED APPLICATIONS	WB, IHC, IF, ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.