

HMP19 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15731	SYNONYMS
SIZE	100 ug	NSG2; Neuron-specific protein family member 2; Protein p19; Hmp19
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	NSG2
PROTEIN NAME	Neuron-specific protein family member 2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human HMP19. at AA range: 10-90

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.