

MRP-L48 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15732	SYNONYMS
SIZE	100 ug	MRPL48; CGI-118; HSPC290; 39S ribosomal protein
HOST	Rabbit	L48; mitochondrial; L48mt; MRP-L48

PRODUCT INFORMATION

GENE NAME	MRPL48
PROTEIN NAME	39S ribosomal protein L48 mitochondrial
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human MRP-L48. at AA range: 130-210

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.