

GAAP Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15733	SYNONYMS
SIZE	100 ug	TMBIM4; GAAP; LFG4; CGI-119; Protein lifeguard 4;
HOST	Rabbit	Golgi anti-apoptotic protein; Protein S1R; Transmembrane BAX inhibitor motif-containing protein 4; Z-protein

PRODUCT INFORMATION

GENE NAME	TMBIM4
PROTEIN NAME	Protein lifeguard 4
IMMUNOGEN	Synthesized peptide derived from the Internal region of human GAAP. at AA range: 90-170

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.