

MRP-S33 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15734	SYNONYMS
SIZE	100 ug	MRPS33; CGI-139; PTD003; 28S ribosomal protein
HOST	Rabbit	S33; mitochondrial; MRP-S33; S33mt

PRODUCT INFORMATION

GENE NAME	MRPS33
PROTEIN NAME	28S ribosomal protein S33 mitochondrial
IMMUNOGEN	Synthesized peptide derived from the Internal region of human MRP-S33. at AA range: 30-110

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.