

CDK5 Activator-binding C42 Polyclonal Antibody

OVERVIEW

| | | |
|--------------------|----------|---|
| CATALOG NO. | PTX15735 | SYNONYMS |
| SIZE | 100 ug | CDK5RAP1; C20orf34; CGI-05; HSPC167; CDK5 |
| HOST | Rabbit | regulatory subunit-associated protein 1; CDK5 activator-binding protein C42 |

PRODUCT INFORMATION

| | |
|---------------------|--|
| GENE NAME | CDK5RAP1 |
| PROTEIN NAME | CDK5 regulatory subunit-associated protein 1 |
| IMMUNOGEN | Synthesized peptide derived from the Internal region of human CDK5 Activator-binding C42. at AA range: 220-300 |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|--|
| SPECIES REACTIVITY | Human |
| TESTED APPLICATIONS | IHC,ELISA |
| RECOMMENDED DILUTION | Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications. |
| DELIVERY CONDITIONS | Blue Ice |
| STABILITY & STORAGE | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles |

For research use only.