

ATP5G2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15736	SYNONYMS
SIZE	100 ug	ATP5G2; PSEC0033; ATP synthase lipid-binding
HOST	Rabbit	protein; mitochondrial; ATP synthase proteolipid P2; ATPase protein 9; ATPase subunit c

PRODUCT INFORMATION

GENE NAME	ATP5G2
PROTEIN NAME	ATP synthase lipid-binding protein mitochondrial
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human ATP5G2. at AA range: 10-90

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.