

Endomucin Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15741	SYNONYMS
SIZE	100 ug	EMCN; EMCN2; MUC14; Endomucin; Endomucin-2;
HOST	Rabbit	Gastric cancer antigen Ga34; Mucin-14; MUC-14

PRODUCT INFORMATION

GENE NAME	EMCN
PROTEIN NAME	Endomucin
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Endomucin. at AA range: 180-260

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	IF,ELISA
RECOMMENDED DILUTION	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.