

CYB5R1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15742	SYNONYMS
SIZE	100 ug	CYB5R1; NQO3A2; NADH-cytochrome b5 reductase 1;
HOST	Rabbit	b5R.1; Humb5R2; NAD(P)H:quinone oxidoreductase type 3 polypeptide A2

PRODUCT INFORMATION

GENE NAME	CYB5R1
PROTEIN NAME	NADH-cytochrome b5 reductase 1
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human CYB5R1. at AA range: 20-100

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.