

PR Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15764	SYNONYMS
SIZE	100 ug	PGR; NR3C3; Progesterone receptor; PR; Nuclear
HOST	Rabbit	receptor subfamily 3 group C member 3

PRODUCT INFORMATION

GENE NAME	PGR
PROTEIN NAME	Progesterone receptor
IMMUNOGEN	Synthesized peptide derived from human PR around the non-phosphorylation site of S190.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.