

Acetyl-Histone H2B (K12) Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15821	SYNONYMS
SIZE	100 ug	H2BFS; Histone H2B type F-S; Histone H2B.s; H2B/s
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	H2BFS
PROTEIN NAME	Histone H2B type F-S
IMMUNOGEN	Synthesized peptide derived from human Histone H2B around the acetylation site of K12. at AA range: 10-90

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Monkey
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.