

Nop132 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15875	SYNONYMS
SIZE	100 ug	NOL8; C9orf34; NOP132; Nucleolar protein 8;
HOST	Rabbit	Nucleolar protein Nop132

PRODUCT INFORMATION

GENE NAME	NOL8
PROTEIN NAME	Nucleolar protein 8
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Nop132. at AA range: 1110-1190

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.