

## PTG Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX15885	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	PPP1R3C; PPP1R5; Protein phosphatase 1 regulatory subunit 3C; Protein phosphatase 1 regulatory subunit 5;
<b>HOST</b>	Rabbit	PP1 subunit R5; Protein targeting to glycogen; PTG

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	PPP1R3C
<b>PROTEIN NAME</b>	Protein phosphatase 1 regulatory subunit 3C
<b>IMMUNOGEN</b>	Synthesized peptide derived from the Internal region of human PTG. at AA range: 20-100

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human,Monkey
<b>TESTED APPLICATIONS</b>	WB,IHC,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**