

LPAAT- ϵ Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15904	SYNONYMS
SIZE	100 ug	AGPAT5; 1-acyl-sn-glycerol-3-phosphate
HOST	Rabbit	acyltransferase epsilon; 1-acylglycerol-3-phosphate O-acyltransferase 5; 1-AGP acyltransferase 5; 1-AGPAT 5; Lysophosphatidic acid acyltransferase epsilon; LPAAT-epsilon

PRODUCT INFORMATION

GENE NAME	AGPAT5
PROTEIN NAME	1-acyl-sn-glycerol-3-phosphate acyltransferase epsilon
IMMUNOGEN	Synthesized peptide derived from the Internal region of human LPAAT- ϵ . at AA range: 210-290

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.