

PLC-XD1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15905	SYNONYMS
SIZE	100 ug	PLCXD1; PI-PLC X domain-containing protein 1
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	PLCXD1
PROTEIN NAME	PI-PLC X domain-containing protein 1
IMMUNOGEN	Synthesized peptide derived from the Internal region of human PLC-XD1. at AA range: 160-240

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.