

## CD172g Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX15909	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	SIRPG; SIRPB2; Signal-regulatory protein gamma;
<b>HOST</b>	Rabbit	SIRP-gamma; CD172 antigen-like family member B; Signal-regulatory protein beta-2; SIRP-b2; SIRP-beta-2; CD antigen CD172g

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	SIRPG
<b>PROTEIN NAME</b>	Signal-regulatory protein gamma
<b>IMMUNOGEN</b>	Synthesized peptide derived from the Internal region of human CD172g. at AA range: 70-150

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human, Monkey
<b>TESTED APPLICATIONS</b>	WB,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**