

ARF GAP1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15935	SYNONYMS
SIZE	100 ug	ARFGAP1; ARF1GAP; ADP-ribosylation factor
HOST	Rabbit	GTPase-activating protein 1; ARF GAP 1; ADP-ribosylation factor 1 GTPase-activating protein; ARF1 GAP; ARF1-directed GTPase-activating protein

PRODUCT INFORMATION

GENE NAME	ARFGAP1
PROTEIN NAME	ADP-ribosylation factor GTPase-activating protein 1
IMMUNOGEN	Synthesized peptide derived from the Internal region of human ARF GAP1. at AA range: 140-220

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.