

GPRC5C Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX15981	SYNONYMS
SIZE	100 ug	GPRC5C; RAIG3; PSEC0087; G-protein coupled
HOST	Rabbit	receptor family C group 5 member C; Retinoic acid-induced gene 3 protein; RAIG-3

PRODUCT INFORMATION

GENE NAME	GPRC5C
PROTEIN NAME	G-protein coupled receptor family C group 5 member C
IMMUNOGEN	Synthesized peptide derived from the Internal region of human GPRC5C. at AA range: 30-110

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.