

Bad Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16061	SYNONYMS
SIZE	100 ug	BAD; BBC6; BCL2L8; Bcl2 antagonist of cell death;
HOST	Rabbit	BAD; Bcl-2-binding component 6; Bcl-2-like protein 8; Bcl2-L-8; Bcl-XL/Bcl-2-associated death promoter

PRODUCT INFORMATION

GENE NAME	BAD
PROTEIN NAME	Bcl2 antagonist of cell death
IMMUNOGEN	Synthesized peptide derived from human Bad around the non-phosphorylation site of S134.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.