

## MYL7 Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX16143	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	MYL7; MYL2A; MYLC2A; Myosin regulatory light chain 2; atrial isoform; MLC-2a; MLC2a; Myosin light chain 2a; Myosin regulatory light chain 7
<b>HOST</b>	Rabbit	

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	MYL7
<b>PROTEIN NAME</b>	Myosin regulatory light chain 2 atrial isoform
<b>IMMUNOGEN</b>	Synthesized peptide derived from the Internal region of human MYL7. at AA range: 80-160

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human,Mouse,Rat,Monkey
<b>TESTED APPLICATIONS</b>	WB,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**