

## Bcl-6 Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX16232	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	BCL6; BCL5; LAZ3; ZBTB27; ZNF51; B-cell
<b>HOST</b>	Rabbit	lymphoma 6 protein; BCL-6; B-cell lymphoma 5 protein; BCL-5; Protein LAZ-3; Zinc finger and BTB domain- containing protein 27; Zinc finger protein 51

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	BCL6
<b>PROTEIN NAME</b>	B-cell lymphoma 6 protein
<b>IMMUNOGEN</b>	Synthesized peptide derived from human Bcl-6 around the non-phosphorylation site of S333.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human,Mouse,Rat
<b>TESTED APPLICATIONS</b>	WB,IHC,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**