

Ribosomal Protein L37 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16254	SYNONYMS
SIZE	100 ug	RPL37; 60S ribosomal protein L37; G1.16
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	RPL37
PROTEIN NAME	60S ribosomal protein L37
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Ribosomal Protein L37. at AA range: 40-120

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.