

Calgizzarin Polyclonal Antibody

OVERVIEW

| | | |
|--------------------|----------|--------------------------------------------------------------------------------------------------------------------|
| CATALOG NO. | PTX16297 | SYNONYMS |
| SIZE | 100 ug | S100A11; MLN70; S100C; Protein S100-A11; |
| HOST | Rabbit | Calgizzarin; Metastatic lymph node gene 70 protein; MLN 70; Protein S100-C; S100 calcium-binding protein A11 |

PRODUCT INFORMATION

| | |
|---------------------|---------------------------------------------------------------------------------------------------|
| GENE NAME | S100A11 |
| PROTEIN NAME | Protein S100-A11 |
| IMMUNOGEN | Synthesized peptide derived from the Internal region of human Calgizzarin. at AA range: 30-110 |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| SPECIES REACTIVITY | Human,Mouse,Rat |
| TESTED APPLICATIONS | IHC,ELISA |
| RECOMMENDED DILUTION | Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications. |
| DELIVERY CONDITIONS | Blue Ice |
| STABILITY & STORAGE | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles |

For research use only.