

Syndecan-2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16317	SYNONYMS
SIZE	100 ug	SDC2; HSPG1; Syndecan-2; SYND2; Fibroglycan;
HOST	Rabbit	Heparan sulfate proteoglycan core protein; HSPG; CD antigen CD362

PRODUCT INFORMATION

GENE NAME	SDC2
PROTEIN NAME	Syndecan-2
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Syndecan-2. at AA range: 50-130

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.