

hCAP-G Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16335	SYNONYMS
SIZE	100 ug	NCAPG; CAPG; NYMEL3; Condensin complex subunit
HOST	Rabbit	3; Chromosome-associated protein G; Condensin subunit CAP-G; hCAP-G; Melanoma antigen NY-MEL-3; Non- SMC condensin I complex subunit G; XCAP-G homolog

PRODUCT INFORMATION

GENE NAME	NCAPG
PROTEIN NAME	Condensin complex subunit 3
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human hCAP-G. at AA range: 920-1000

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.