

Cerebellin 3 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16350	SYNONYMS
SIZE	100 ug	CBLN3; Cerebellin-3
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	CBLN3
PROTEIN NAME	Cerebellin-3
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Cerebellin 3. at AA range: 100-180

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.