

Spectrin β II Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16444	SYNONYMS
SIZE	100 ug	SPTBN1; SPTB2; Spectrin beta chain; non-erythrocytic
HOST	Rabbit	1; Beta-II spectrin; Fodrin beta chain; Spectrin, non-erythroid beta chain 1

PRODUCT INFORMATION

GENE NAME	SPTBN1
PROTEIN NAME	Spectrin beta chain non-erythrocytic 1
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Spectrin β II. at AA range: 620-700

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.