

SSTR1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16473	SYNONYMS
SIZE	100 ug	SSTR1; Somatostatin receptor type 1; SS-1-R; SS1-R;
HOST	Rabbit	SS1R; SRIF-2

PRODUCT INFORMATION

GENE NAME	SSTR1
PROTEIN NAME	Somatostatin receptor type 1
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human SSTR1. at AA range: 20-100

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.