

TAF II p18 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16518	SYNONYMS
SIZE	100 ug	TAF13; TAF2K; TAFII18; Transcription initiation factor
HOST	Rabbit	TFIID subunit 13; Transcription initiation factor TFIID 18 kDa subunit; TAF(II)18; TAFII-18; TAFII18

PRODUCT INFORMATION

GENE NAME	TAF13
PROTEIN NAME	Transcription initiation factor TFIID subunit 13
IMMUNOGEN	Synthesized peptide derived from the Internal region of human TAF II p18. at AA range: 40-120

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.