

TAT Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16526	SYNONYMS
SIZE	100 ug	TAT; Tyrosine aminotransferase; TAT; L-tyrosine:2-oxoglutarate aminotransferase
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	TAT
PROTEIN NAME	Tyrosine aminotransferase
IMMUNOGEN	Synthesized peptide derived from the Internal region of human TAT.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	IHC,ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.