

TR α Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16571	SYNONYMS
SIZE	100 ug	THRA; EAR7; ERBA1; NR1A1; THRA1; THRA2;
HOST	Rabbit	Thyroid hormone receptor alpha; Nuclear receptor subfamily 1 group A member 1; V-erbA-related protein 7; EAR-7; c-erbA-1; c-erbA-alpha

PRODUCT INFORMATION

GENE NAME	THRA
PROTEIN NAME	Thyroid hormone receptor alpha
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human TR α . at AA range: 10-90

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	WB,IHC,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.