

p53 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16616	SYNONYMS
SIZE	100 ug	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-
HOST	Rabbit	13; Phosphoprotein p53; Tumor suppressor p53

PRODUCT INFORMATION

GENE NAME	TP53
PROTEIN NAME	Cellular tumor antigen p53
IMMUNOGEN	Synthesized peptide derived from human p53 around the non-acetylation site of K317. at AA range: 260-340

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.