

TSSC1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16636	SYNONYMS
SIZE	100 ug	TSSC1; Protein TSSC1; Tumor-suppressing STF cDNA
HOST	Rabbit	1 protein; Tumor-suppressing subchromosomal transferable fragment candidate gene 1 protein

PRODUCT INFORMATION

GENE NAME	TSSC1
PROTEIN NAME	Protein TSSC1
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human TSSC1. at AA range: 260-340

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.