

THP Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16659	SYNONYMS
SIZE	100 ug	UMOD; Uromodulin; Tamm-Horsfall urinary
HOST	Rabbit	glycoprotein; THP

PRODUCT INFORMATION

GENE NAME	UMOD
PROTEIN NAME	Uromodulin
IMMUNOGEN	Synthesized peptide derived from the Internal region of human THP.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.