

ZnT-4 Polyclonal Antibody

OVERVIEW

| | | |
|--------------------|----------|---|
| CATALOG NO. | PTX16734 | SYNONYMS |
| SIZE | 100 ug | SLC30A4; ZNT4; Zinc transporter 4; ZnT-4; Solute carrier family 30 member 4 |
| HOST | Rabbit | |

PRODUCT INFORMATION

| | |
|---------------------|--|
| GENE NAME | SLC30A4 |
| PROTEIN NAME | Zinc transporter 4 |
| IMMUNOGEN | Synthesized peptide derived from the Internal region of human ZnT-4. at AA range: 40-120 |

APPLICATIONS AND STORAGE

| | |
|--------------------------------|--|
| SPECIES REACTIVITY | Human,Mouse,Rat |
| TESTED APPLICATIONS | WB,IHC,ELISA |
| RECOMMENDED DILUTION | Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications. |
| DELIVERY CONDITIONS | Blue Ice |
| STABILITY & STORAGE | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles |

For research use only.