

SgK269 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16778	SYNONYMS
SIZE	100 ug	PEAK1; KIAA2002; Pseudopodium-enriched atypical
HOST	Rabbit	kinase 1; Sugen kinase 269; Tyrosine-protein kinase SgK269

PRODUCT INFORMATION

GENE NAME	PEAK1
PROTEIN NAME	Pseudopodium-enriched atypical kinase 1
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human SgK269. at AA range: 80-160

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.