

CSRNP3 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16791	SYNONYMS
SIZE	100 ug	CSRNP3; FAM130A2; TAIP2; Cysteine/serine-rich
HOST	Rabbit	nuclear protein 3; CSRNP-3; Protein FAM130A2; TGF-beta-induced apoptosis protein 2; TAIP-2

PRODUCT INFORMATION

GENE NAME	CSRNP3
PROTEIN NAME	Cysteine/serine-rich nuclear protein 3
IMMUNOGEN	Synthesized peptide derived from the Internal region of human CSRNP3. at AA range: 140-220

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.