

GPR68 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16813	SYNONYMS
SIZE	100 ug	GPR68; OGR1; Ovarian cancer G-protein coupled receptor 1; OGR-1; G-protein coupled receptor 68; GPR12A; Sphingosylphosphorylcholine receptor
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	GPR68
PROTEIN NAME	Ovarian cancer G-protein coupled receptor 1
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human GPR68. at AA range: 250-330

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse
TESTED APPLICATIONS	WB,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.