

## Olfactory receptor 2C3 Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX16821	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	OR2C3; OR2C4; OR2C5P; Olfactory receptor 2C3;
<b>HOST</b>	Rabbit	Olfactory receptor 2C4; Olfactory receptor 2C5; Olfactory receptor OR1-30

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	OR2C3
<b>PROTEIN NAME</b>	Olfactory receptor 2C3
<b>IMMUNOGEN</b>	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2C3. at AA range: 200-280

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human
<b>TESTED APPLICATIONS</b>	WB,IF,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**