

CaMKII $\beta/\gamma/\delta$ Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16828	SYNONYMS
SIZE	100 ug	CAMK2B; CAM2; CAMK2; CAMKB;
HOST	Rabbit	Calcium/calmodulin-dependent protein kinase type II subunit beta; CaM kinase II subunit beta; CaMK-II subunit beta; CAMK2G; CAMK; CAMK-II; CAMKG; Calcium/calmodulin-dependent protein kinase type II subunit gamma;

PRODUCT INFORMATION

GENE NAME	CAMK2B
PROTEIN NAME	Calcium/calmodulin-dependent protein kinase type II subunit beta
IMMUNOGEN	Synthesized peptide derived from human CaMKII $\beta/\gamma/\delta$ around the non-phosphorylation site of T287.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse, Rat
TESTED APPLICATIONS	WB, IHC, ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.