

Frizzled-4 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16848	SYNONYMS
SIZE	100 ug	FZD4; Frizzled-4; Fz-4; hFz4; FzE4; CD antigen CD344
HOST	Rabbit	

PRODUCT INFORMATION

GENE NAME	FZD4
PROTEIN NAME	Frizzled-4
IMMUNOGEN	Synthesized peptide derived from the Internal region of human Frizzled-4. at AA range: 100-180

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Mouse,Rat
TESTED APPLICATIONS	IF,ELISA
RECOMMENDED DILUTION	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.