

Caspase-1 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16859	SYNONYMS
SIZE	100 ug	CASP1; IL1BC; IL1BCE; Caspase-1; CASP-1;
HOST	Rabbit	Interleukin-1 beta convertase; IL-1BC; Interleukin-1 beta-converting enzyme; ICE; IL-1 beta-converting enzyme; p45

PRODUCT INFORMATION

GENE NAME	CASP1
PROTEIN NAME	Caspase-1
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Caspase-1. at AA range: 270-350

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.