

## MAD1 Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX16894	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	MAD1L1; MAD1; TXBP181; Mitotic spindle assembly
<b>HOST</b>	Rabbit	checkpoint protein MAD1; Mitotic arrest deficient 1-like protein 1; MAD1-like protein 1; Mitotic checkpoint MAD1 protein homolog; HsMAD1; hMAD1; Tax- binding protein 181

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	MAD1L1
<b>PROTEIN NAME</b>	Mitotic spindle assembly checkpoint protein MAD1
<b>IMMUNOGEN</b>	Synthesized peptide derived from human MAD1 around the non-phosphorylation site of S428.

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human
<b>TESTED APPLICATIONS</b>	IF,ELISA
<b>RECOMMENDED DILUTION</b>	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**