

Nck-2 Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16938	SYNONYMS
SIZE	100 ug	NCK2; GRB4; Cytoplasmic protein NCK2; Growth
HOST	Rabbit	factor receptor-bound protein 4; NCK adaptor protein 2; Nck-2; SH2/SH3 adaptor protein NCK-beta

PRODUCT INFORMATION

GENE NAME	NCK2
PROTEIN NAME	Cytoplasmic protein NCK2
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human Nck-2. at AA range: 300-380

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human, Mouse
TESTED APPLICATIONS	IHC, ELISA
RECOMMENDED DILUTION	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.