

## NFATc2IP Polyclonal Antibody

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX16967	<b>SYNONYMS</b>
<b>SIZE</b>	100 ug	NFATC2IP; NIP45; NFATC2-interacting protein; 45
<b>HOST</b>	Rabbit	kDa NF-AT-interacting protein; 45 kDa NFAT-interacting protein; Nuclear factor of activated T-cells; cytoplasmic 2-interacting protein

### PRODUCT INFORMATION

---

<b>GENE NAME</b>	NFATC2IP
<b>PROTEIN NAME</b>	NFATC2-interacting protein
<b>IMMUNOGEN</b>	Synthesized peptide derived from the Internal region of human NFATc2IP. at AA range: 120-200

### APPLICATIONS AND STORAGE

---

<b>SPECIES REACTIVITY</b>	Human
<b>TESTED APPLICATIONS</b>	WB,ELISA
<b>RECOMMENDED DILUTION</b>	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
<b>DELIVERY CONDITIONS</b>	Blue Ice
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

**For research use only.**