

MASTL Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16969	SYNONYMS
SIZE	100 ug	MASTL; GW; GWL; THC2; Serine/threonine-protein
HOST	Rabbit	kinase greatwall; GW; GWL; hGWL; Microtubule-associated serine/threonine-protein kinase-like; MAST-L

PRODUCT INFORMATION

GENE NAME	MASTL
PROTEIN NAME	Serine/threonine-protein kinase greatwall
IMMUNOGEN	Synthesized peptide derived from the C-terminal region of human MASTL.

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human
TESTED APPLICATIONS	WB,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.