

Atg4C Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX16971	SYNONYMS
SIZE	100 ug	ATG4C; APG4C; AUTL1; AUTL3; Cysteine protease
HOST	Rabbit	ATG4C; AUT-like 3 cysteine endopeptidase; Autophagin-3; Autophagy-related cysteine endopeptidase 3; Autophagy-related protein 4 homolog C

PRODUCT INFORMATION

GENE NAME	ATG4C
PROTEIN NAME	Cysteine protease ATG4C
IMMUNOGEN	Synthesized peptide derived from the N-terminal region of human Atg4C. at AA range: 40-120

APPLICATIONS AND STORAGE

SPECIES REACTIVITY	Human,Monkey
TESTED APPLICATIONS	WB,IHC,IF,ELISA
RECOMMENDED DILUTION	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
DELIVERY CONDITIONS	Blue Ice
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles

For research use only.